

Liposuction Cannulas

One Piece Cannula special
\$145.00 ea.



AMT offers the highest quality stainless steel cannulas at a competitive price point. Each cannula is precisely machined, polished and assembled to your exact specifications.

We have included some of the more common cannula choices in our catalog. However, *Please let us know if there is something that you would like that you don't see here.*

To Order a Cannula, please specify the following:

1. **Tip Shape**
2. **Tip Style**
3. **Diameter**
4. **Length**
5. **Handle Type**
6. **Straight or Curved**
7. **Vent Hole on Handle** Yes or No?
8. **Perfect Alignment** to indicate orientation of cannula. I.E. For one sided orifice opening.
9. **Autoclave Tray or Set** Please let us know if you would like to put together an entire set?

Cannula Tips



Bullet Tip



Blunt Tip



Spatula Tip



Wedge Tip

Cannula Styles



Mercedes



Double Mercedes



Accelerator



Double Accelerator



Illouz Standard



2-Hole Standard



3-Hole Standard



1-Hole Spatula



Texas



1-Hole V-Notched



Cobra K



Cobra



Cobra +




Cobra K +



Capri



Grater-Round



Sattler-Blugerman (24 Holes, 4 rows)



Sattler-Blugerman Draining Cannula (12 Holes)

Liposuction Cannulas cont.

Diameters Outside Diameters (not to scale)

- 1.2mm ◦ 18g
- 1.4mm ◦ 17g
- 1.6mm ◦ 16g
- 2mm ◦ 14g
- 2.7mm ◦ 12g
- 3mm ◦ 11g
- 3.4mm ◦ 10g
- 4mm ◦ 8g
- 5mm ◦ 6g

Standard Lengths

7cm - 9cm - 10cm - 15cm - 23cm - 30cm - 35cm

Curve Preference Choice



Straight - Convex Curve - Concave - Swan Neck (Pictured here)

Cannula Handle Options

We have included some of the more common cannula handle choices in our catalog. *However, please let us know if there is something that you would like that you don't see here.*

One Piece Aluminum Handles, Contour Handles (With half the weight and more control), Sure Fit Luer lock handles, Two Piece Quick Screw (Delrin) handles and NEW One Piece Handles with a swivel barb to allow easier turning of the handle. (A Tube Swivel Adaptor is also pictured)



NEW One-Piece Swivel Barb Handle



NEW Tube Swivel Adaptor

Thank you for your interest in our Cannula product offering. Please contact an AMT Sales representative to help with a custom order. Again, please let us know if there is something that you would like that you don't see here.