

Tumescent Liposuction

Specials


Infiltration/Infusion Needles with Sure Fit Luer Lock

\$80 ea.

Available in single, multi-hole and Garden hose style versions. Infusion cannulas are strong, slightly flexible and resilient (all on Luer Lock Hubs). Standard diameters are 1.5mm, 2.0mm, and 3.0mm. Standard lengths are 17cm, 25cm and 32cm.

Custom Configurations are available.

Standard Spiral (6 Hole) & Garden Hose Style (4 Holes at 90 Degrees) Shown Here.



Liposuction Infusion Handle Options:



INFUSION HANDLE WITH ON/OFF CONTROL LUER LOCK CONNECT

\$ 175 ea.

The on/off valve workings have been improved by sealing systems that will better withstand the high autoclave temperatures. This improved handle allows for fingertip control during the infusion process. Lightweight Aluminum one piece design.

Order No. **AMT1040-INFU-C**



INFUSION HANDLE WITH STOP COCK CONTROL

\$ 160 ea.

The level handle shows flow direction for easy visual checks, quarter turn shot-off. The valve is made of brass and chrome-plated finish. The stop cock control can be rotated to any direction by releasing the handle cap and tight the cap back after you choose what direction you want. Two Piece design. May be used with out stopcock.

Order No. **AMT1540-INFU-SCC**