

AMT Body Liposuction SET

Designed for Body aspiration procedures

Tip Styles & Sizes

- (1) 3-Port Radial (Mercedes)
- (1) 3-Port Radial (Mercedes) 3mm x 23cm
- (1) 3-Port Radial (Mercedes) 4mm x 30cm
- (1) Trip Port- Accelerator 3mm x 23cm
- (1) Tri Port -Accelerator 4mm x 30cm



Handle Options:

- A) All One piece handles
- B) Two Piece - Luer Lock Handles - Qty. 2
- C) Two Piece - Delrin Handles - Qty. 2


Infusion Needles:

- (1) Infusion Handle, luer lock, with on/off control
- (1) 3.0mm x 23cm spiral 6 hole fusion needle, luer lock
- (1) 2.0mm x 17cm spiral 6 hole infusion needle, luer lock



Item No. AMT-BodyLipo

- (1) Autoclavable Tray Optional



Helical & K-Wedge - Body Liposuction Set

Helical Pattern Single Slot Cannulas

- 16g X 8cm
- 16g X 12cm
- 14g X 15cm
- 14g X 23cm
- 12g X 15cm
- 12g X 9cm


K-Wedge Cannulas

- 14g X 10cm
- 14g X 15cm
- 12g X 10cm
- 12g X 15cm

Autoclavable Tray



Helical Pattern Single Slot compare to **Capistrano™** style. Multiple round holes in helical pattern with one single slot distally.



K-Wedge compare to **Finesse™** Style. Two holes on one side and a spatula-like tip. Often used near the dermis during last stages Of liposuction.