




800-268-2937

Custom Design Cannulas



Standard One Piece  
**\$100.00**



Swivel Barb One Piece  
**\$115.00**



Delrin Screw Style  
**\$95.00**



Leur Lock Cannulas  
**\$85.00**

**Custom Design Steps:** Please Indicate your preferences for the following:  
We will machine, polish & assemble to your specifications.

**1.) Tip Shape**



Bullet Tip



Blunt Tip



Spatula Tip



Wedge Tip

**2.) Tip Style** - only limited by available surface area and your imagina-

**(Click Tip Style Category or Scroll through PDF for Popular Selections)**



**3.) Standard Outside Diameters**

MicroCannulas

1.2mm (18g); 1.4mm (17g); 1.6mm (16g); 2mm (14g); 2.7mm (12g)

Body Cannulas

3mm (11g); 3.4mm (10g); 4mm (8g); 5mm (6g)

**4.) Standard Lengths** 7cm - 9cm - 10cm - 15cm - 23cm - 30cm -

**5.) Curve Preference**



Straight



Convex Curve



Swan Neck or S-Curve



Concave Curve

## Conservative Cannula Tips

[Back to Top](#) ▶

**Popular Item**

One Hole Standard Single Slot is located 180 degrees from thumb rest on the hub.



Two Hole Standard Often used in face and neck for sculpting and helping to provoke skin retraction.



Three Hole Standard



**Popular Item**

Mercedes circumferential three port pattern near the distal tip of shaft. Multi-Plane soft tissue harvesting.



**Popular Item**

Triport compare to **Accelerator "III"™** style. Three holes in a triangular pattern.



**Popular Item**

Pyramid compare to **Las Vegas™** style A single distal opening with two proximal openings. Ideal for feathering.



**Popular Item**

Spiral Cannula compare to **Candy Cane™** style. Features three elongated openings in a spiral pattern.



## Spatula & Wedge Tips

[Back to Top](#) ▶

**Popular Item**

One Hole Spatula Named for its flattened profile, the spatula is ideal for cervical and facial procedures.



**Popular Item**

Helical Pattern Single Slot compare to **Capistrano™** style. Multiple round holes in helical pattern with one single slot distally.



K-Wedge compare to **Finesse™** Style. Two holes on one side and a spatula-like tip. Often used near the dermis during last stages of liposuction.



## Draining Cannula

[Back to Top](#) ▶

Draining Cannula (12 holes)



## Superficial Cannula Tips

[Back to Top](#) ▶

Bevel tips works well for needle extraction of soft tissue and as a soft tissue reinjection needle.



One Hole "V" Notched Compare to **Toledo V™**.



Etching Cannula Compare to **Gilliland™** 4 openings in same plane. Distal ends of openings are raised for etching of subcutaneous tissue.



Texas



Texas+




## Dissecting Cannula Tips

[Back to Top](#) ▶

Flare Mercedes Compare to **Becker Basket™**  
Riblike projections facilitate breakdown of tissue prior to aspiration, increasing speed and efficiency in a wide range of procedures.



24 hole 4 Rows) Style I, II or III  
Compare to **Sattler™** or **Blugerman** style.



Grater Round Compare to **Becker Grater Round™**



3 Hole V Notch Dissector Compare to **Aiache™**



CRD - Cellulite Release Dissector



## Aggressive Cannula Tips

[Back to Top](#) ▶

Cobra



Cobra +



Cobra K



Cobra K +



Keel Cobra

