

AMT Face & Neck Liposuction SET

Designed for small area aspiration procedures

Tip Styles & Sizes

- (1) Accelerator 2mm x 15cm
- (1) Accelerator 2mm x 23cm
- (1) Accelerator 3mm x 15cm
- (1) Spatula 2mm x 15cm
- (1) Spatula 3mm x 15cm



Handle Options:

- A) All One piece handles
- B) Two Piece - Luer Lock Handles - Qty. 2
- C) Two Piece - Delrin Handles - Qty. 2


Infusion Needles:

- (1) Infusion Handle, luer lock, with on/off control
- (1) 1.5mm x 15cm spatula infusion needle, luer lock
- (1) 1.5mm x 10cm spatula infusion needle, luer lock



Item No. AMT-FaceLipo

- (1) Autoclavable Tray Optional



Helical & K-Wedge - Facial Liposuction Set

Helical Pattern Single Slot Cannulas

20g X 6cm

18g X 5cm

16g X 8cm

K-Wedge Cannulas

16g X 8cm

14g X 10cm

Autoclavable Tray



Helical Pattern Single Slot compare to **Capistrano™** style. Multiple round holes in helical pattern with one single slot distally.



K-Wedge compare to **Finesse™** Style. Two holes on one side and a spatula-like tip. Often used near the dermis during last stages Of liposuction.

