



Helical Pattern Single Slot compare to **Capistrano™** style. Designed to maximize liposuction efficiency, openings are arranged along entire cannula surface, some are always directed toward the skin. Multiple round holes in helical pattern with one single slot distally. Distal Single slot used to facilitate cleaning and removal of debris.



**\$85.00** each

K-Wedge compare to **Finesse™** Style. Designed to minimize risk of injury to the deep surface of the dermis when liposuction is done near the skin. Two holes on one side and a spatula-like tip. Often used during last stages of liposuction when dermis is more vulnerable to trauma.



Standard Diameters

MicroCannulas

1.2mm (18g); 1.4mm (17g); 1.6mm (16g); 2mm (14g); 2.7mm (12g)

Body Cannulas

3mm (11g); 3.4mm (10g); 4mm (8g); 5mm (6g)

Standard Lengths

7cm, 9cm, 10cm, 15cm, 23cm, 30cm, 35cm

**Helical & K-Wedge** Body Liposuction Set

Helical Pattern Single Slot Cannulas

16g X 8cm

16g X 12cm

14g X 15cm

14g X 23cm

12g X 15cm

12g X 9cm

K-Wedge Cannulas

14g X 10cm

14g X 15cm

12g X 10cm

12g X 15cm

Autoclavable Tray

**Helical & K-Wedge** Facial Liposuction Set

Helical Pattern Single Slot Cannulas

20g X 6cm

18g X 5cm

16g X 8cm

K-Wedge Cannulas

16g X 8cm

14g X 10cm

Autoclavable Tray

