



NEW Laser Assisted Liposuction Instruments

Compatible with all laser liposuction systems including **Smartlipo™**, **Smartlipo MPX™**, **CoolLipo™**, **SlimLipo™**, **ProLipo PLUS™** and others

Laser Fiber Guide Handles & Sheaths are used for guiding the laser fiber to the areas that you would like to melt fat.

Compare & Save

Laser Fiber Guide Handle


Designed to be used with all laser assisted liposuction systems. The Laser Fiber Handle is 5 1/4" long & is available in 3 diameters. 1/2", 3/4" & 1".



Laser Fiber Guide Handles
\$195.00

Laser Fiber Sheaths

Luer lock Sheaths to be used with laser assisted liposuction systems. Cannulas come in a variety of sizes to accommodate any systems' fiber guide.




Laser Fiber Sheaths
\$95.00

NEW Laser Assist Liposuction Handle

Allows for simultaneous melting & suction of fat, for greater control and accuracy of fat removal. Provides user with immediate feedback to help determine endpoints & decrease procedure time.

"Our Study of 50 cases showed, the Guide gave superior results over the traditional technique (separate lasing and suction) while significantly reducing procedure times"

Bruce E. Katz, MD
New York, NY



Combo Laser/Lipo Handle

Standard Diameter 2mm Facial ; 3.4mm Body

Standard Lengths 10cm, 15cm, 20cm or 25cm

Laser Assist Lipo Handle
Contact for Quote