




Luer Lock Handle Options

Luer Lock Handle Lightweight Aluminum. Designed for all standard luer lock aspiration cannulas. Large 1" or Small 3/4" Diam.
\$85.00ea



Luer Lock Handle w/ Swivel Barb Designed for all standard luer lock aspiration cannulas. Enables easier turning of handle during aspiration. 3/4" Diam
\$115.00ea



- 1.) Design Your Own Luer Lock Cannula
- Or
- 2.) Choose from our most popular tip styles

We will machine, polish & assemble to your specifications.



\$85.00 ALL Luer Lock Cannula Styles

Tip Style - only limited by available surface area and your imagination.

(Click Tip Style Category or Scroll through PDF for Popular Selections)

Conservative Cannulas ▶	Spatula & Wedge Tips ▶	Draining Cannulas ▶
Dissecting Cannulas ▶	Superficial Cannulas ▶	Aggressive Cannulas ▶

Standard Outside Diameters

MicroCannulas

1.2mm (18g); 1.4mm (17g); 1.6mm (16g); 2mm (14g); 2.7mm (12g)

Body Cannulas

3mm (11g); 3.4mm (10g); 4mm (8g); 5mm (6g)

Standard Lengths 7cm - 9cm - 10cm - 15cm - 23cm - 30cm - 35cm

Curve Preference



Straight



Convex Curve



Swan Neck or S-Curve



Concave Curve

Conservative Cannula Tips

[Back to Top](#) ▶

Popular Item

One Hole Standard Single Slot is located 180 degrees from thumb rest on the hub.



Two Hole Standard Often used in face and neck for sculpting and helping to provoke skin retraction.



Three Hole Standard



Popular Item

Mercedes circumferential three port pattern near the distal tip of shaft. Multi-Plane soft tissue harvesting.



Popular Item

Triport compare to **Accelerator "III"™** style. Three holes in a triangular pattern.



Popular Item

Pyramid compare to **Las Vegas™** style A single distal opening with two proximal openings. Ideal for feathering.



Popular Item

Spiral Cannula compare to **Candy Cane™** style. Features three elongated openings in a spiral pattern.



Spatula & Wedge Tips

[Back to Top](#) ▶

Popular Item

One Hole Spatula Named for its flattened profile, the spatula is ideal for cervical and facial procedures.



Popular Item

Helical Pattern Single Slot compare to **Capistrano™** style. Multiple round holes in helical pattern with one single slot distally.



K-Wedge compare to **Finesse™** Style. Two holes on one side and a spatula-like tip. Often used near the dermis during last stages of liposuction.



Draining Cannula

[Back to Top](#) ▶

Draining Cannula (12 holes)



Superficial Cannula Tips

Bevel tips works well for needle extraction of soft tissue and as a soft tissue reinjection needle.


One Hole "V" Notched Compare to **Toledo V™**.

Etching Cannula Compare to **Gilliland™** 4 openings in same plane. Distal ends of openings are raised for etching of subcutaneous tissue.

Texas

Texas+

[Back to Top](#) ▶



Dissecting Cannula Tips

Flare Mercedes Compare to **Becker Basket™**
Riblike projections facilitate breakdown of tissue prior to aspiration, increasing speed and efficiency in a wide range of procedures.

24 hole 4 Rows) Style I, II or III
Compare to **Sattler™** or **Blugerman** style.

Grater Round Compare to **Becker Grater Round™**

3 Hole V Notch Dissector Compare to **Aiache™**

CRD - Cellulite Release Dissector

[Back to Top](#) ▶



Aggressive Cannula Tips

Cobra

Cobra +

Cobra K

Cobra K +

Keel Cobra

[Back to Top](#) ▶

