



800-268-2937

Micro-Injection Guns

Micro-Injection Gun The new modification gun design for small tissue injection with extra fine teeth. It is three times smaller per squeeze compared to the previous designs. See chart for the approximate result per squeeze when you use it with different sized syringes. This set travel one teeth per squeeze.

B-D Syringe	result per squeeze in mL
3mL	.11 mL
5mL	.20 mL
10mL	.30 mL
60mL	3.21264ml



1, 3, 5 & 10cc Injection Gun



60 cc Injection Gun

Call **800-268-2937** or email info@AestheticMT.com
for a **Customized Quote**.